

**Indian Institute of Technology (IIT) Bombay, Powai
Mumbai-400076**

**Department of Biosciences & Bioengineering
Phone: 022-25764757**

**Registration form for Surface Plasmon Resonance Facility
(spr.bios@iitb.ac.in)**

Date: _____

External Registration Number (Office use only):

Name of User:

**Name of
Institution/Organization:** _____

**Name of the
Dept/Div/Sec:** _____

Email and Tel.No.: _____

**Nature of samples
involved:** _____

**Number of
Samples to be
tested:** _____

Type of Analysis required: _____

**Molecular weight of ligand and analytes
involved:** _____

Isoelectric point (pI) of the ligand involved: _____

Available concentration of ligand and analytes: _____

In case of small molecule analyte: Is the molecule dissolved in organic solvent? Is yes, please specify: _____

Does your analyte tend to aggregate? : _____

Any earlier experimental indications of the ligand-analyte interaction: _____

Conformational changes on ligand-analyte interaction: _____

Any specific temperature requirement: _____

Any other details to be shared: _____

Kindly mention details (and bring along for discussion) from literature search performed on similar studies:

***Ligand: The ligand is the interaction partner attached to the surface.**

***Analyte: The analyte is the interaction partner that passes in solution over the surface.**

INSTRUCTIONS FOR SAMPLE PREPARATION

- Experiments should be discussed with the facility in-charge before appointment.
- Purity of samples is extremely important for generating good data.
- Protein concentrations should be measured accurately before starting the experiment.
- The molecular weight as well as the pI of the proteins should be known before immobilization.
- All buffers should be filtered through 0.22 micron filters and degassed.
- Do not degas buffers containing detergent. Add detergent after degassing.
- For organic solvent containing buffer, filter using organic solvent resistant membrane.
- Cell extracts and nanoparticles can block integrated micro fluidic cartridges and syringes.
- Any query regarding your SPR experiment can be emailed on **spr.bios@iitb.ac.in**
- Appointments will be provided as per que and the user will be informed about the same.
- Kindly perform literature review on similar work and accumulate as much information as possible for good quality data.

Whenever the prepared samples are used in the publications appropriate acknowledgement of usage of IIT Bombay SPR facility must be mentioned. The details should be forwarded to spr.bios@iitb.ac.in

GIVEN MATERIAL IS NOT POISONOUS OR TOXIC IN ANY WAY: _____

We agree to acknowledge the Surface Plasmon Resonance (SPR) Central Facility of IIT BOMBAY in our Publications/Reports/Thesis in which the data is used with due feedback through email.

**Signature of
User:** _____

Sample received (date): _____

Sample analysis completion (date): _____

Signature of concerned Staff-in-charge/TA: _____